



Trolleys, Stands, Carts and Storage

Imaging Solutions takes pride in offering a wide range of Trolleys, Stands, Carts, and Storage products. We manufacture a full selection of RADstore carts and also assemble other products here in Australia. The offering high levels of utility and versatility through inbuilt features which leverage modular reconfigurable product design to meet changing procedure needs and customisation preferences. Our products are designed to meet the unique demands of different healthcare settings and are made of high-quality materials.

your single source supplier™





Trusted provider to medical imaging professionals since 2000

Imaging Solutions is celebrating over 25 years as a trusted provider of market-leading medical imaging equipment and accessories. With an unwavering commitment to innovation, quality, and service, we continue to deliver cutting-edge solutions that empower healthcare professionals worldwide.

25
years

your single source supplier[™]

Table of Contents

Stainless Steel Storage	4
Instrument Trolleys	4
Perspex Dressing Trolley	4
Dressing Trolley	5
Flat Top Trolley	6
IV Poles	7
Endoscopy Carts	8
Endoscopy Travel Carts	8
Endoscopy Procedure Carts	10
Carts	11
Roam Carts	12
Pace Carts	22

Regulatory Notice

Imaging Solutions operates in Australia, New Zealand, the United States, the European Union, and engages with resellers and distributors globally. Product availability in this catalogue may vary by region and is subject to local regulatory classifications and approvals. Certain products may only be available in jurisdictions where they have been registered or authorized by the relevant regulatory authorities. For details on product availability in your region, please contact Imaging Solutions.

Instrument Trolleys

Designed for durability and efficiency, our Instrument Trolley provides a secure and reliable solution for organizing and transporting medical instruments. Crafted from high-quality stainless steel with a fully welded construction, it ensures long-term stability in demanding healthcare settings.

Secure Shelving System

Three perimeter rails per shelf keep instruments securely in place, minimizing the risk of spills or displacement in fast-paced environments.

Smooth Mobility

Fitted with 4 x 100mm swivel castors, this trolley allows effortless movement across various surfaces, ensuring seamless workspace transitions.

Robust Stainless Steel Construction

Its fully welded stainless steel frame provides exceptional strength and durability, making it ideal for high-traffic medical settings.

Unmatched Quality

With its sturdy build, secure storage, and effortless maneuverability, the Instrument Trolley is an essential asset for professionals requiring a dependable, high-performance transport solution.



Part Number	EMY-SP256	EMY-SP256.1	EMY-SP256.2	EMY-SP256-3
Width	500	750	1000	500
Depth	500	500	500	500
Height	900	900	900	900
Shelves	2	2	2	3
Rails	3 on each shelf	3 on each shelf	3 on each shelf	3 on each shelf
Castors	100mm Swivel	100mm Swivel	100mm Swivel	100mm Swivel

Perspex Dressing Trolley

Our Perspex Dressing Trolley offers all the durability and functionality of our standard dressing trolley, with the added benefit of clear perspex drawers for easy content visibility and labeling. Precision ball-bearing castors ensure smooth maneuverability with minimal effort, while full-extension ball-bearing runners provide effortless drawer access without the risk of opening during movement. Designed with safety and infection control in mind, it features full safety folds and smooth, curved handles for a secure grip. Built from high-quality stainless steel, this trolley is backed by a 10-Year Workmanship Guarantee, ensuring long-term reliability in demanding medical environments.

- Dimensions: 500w x 500d x 900h
- Drawers: 4 Clear Perspex
- Shelves: 1 w/ Rails
- Castors: 125mm Swivel

Perspex Dressing Trolley

EMY-SS13VP-P



Dressing Trolley

New Curved Ergonomic Stainless Steel / Welded support frame

Designed to meet the evolving demands of medical environments, our Dressing Trolley combines durability, functionality, and ergonomic excellence. Crafted from high-grade stainless steel with a fully welded construction, it ensures long-term reliability in even the most demanding settings. Precision engineering enhances ease of use, while its sleek, easy-to-clean design supports infection control. Available in multiple sizes and drawer configurations, this trolley delivers the perfect balance of strength and practicality to suit diverse healthcare needs.

Premium Construction

Crafted from high-quality 304 stainless steel, this dressing trolley is built for exceptional durability. Fully welded top and bottom surfaces provide long-lasting sturdiness, meeting the demands of fast-paced medical environments. Its design prioritizes easy cleaning and optimal infection control.

Smooth Functionality

Precision ball-bearing castors ensure effortless maneuverability, while full-extension ball-bearing drawer runners enable smooth access to contents. Designed for security, drawers remain firmly closed during movement.

Part Number	EMY-SS10	EMY-SS10.1	EMY-SS10.2
Width	500	750	1000
Depth	500	500	500
Height	900	900	900
Drawers	1	1	1
Shelves	2 w/ Rails	2 w/ Rails	2 w/ Rails
Castors	100mm Swivel	100mm Swivel	100mm Swivel

Part Number	EMY-SS11	EMY-SS11.1	EMY-SS11.2
Width	500	750	1000
Depth	500	500	500
Height	900	900	900
Drawers	2	2	2
Shelves	2 w/ Rails	2 w/ Rails	2 w/ Rails
Castors	100mm Swivel	125mm Swivel	4 x100mm Swivel

Part Number	EMY-SS12	EMY-SS12.1	EMY-SS12.2
Width	500	750	1000
Depth	500	500	500
Height	900	900	900
Drawers	3	3	3
Shelves	2 w/ Rails	2 w/ Rails	2 w/ Rails
Castors	125mm Swivel	125mm Swivel	125mm Swivel

Safety & Ergonomics

Engineered with medical environments in mind, the lower shelf is set at 300mm from the floor to support occupational health and safety requirements, reducing bending and lifting strain. Full safety folds come standard, ensuring secure and reliable use.

Unmatched Quality

Our dressing trolley stands out with fully seam-welded construction, delivering exceptional strength and longevity. This meticulous design also enhances infection control by ensuring seamless, easy-to-clean surfaces.



Drawer Colour Customisations

Our Dressing Trolleys can be customized to suit specific medical needs, with options such as mounted IV poles for added functionality and drawer color variations for easy identification and organization.



Trolley Mounted IV Pole
EMY-AO-IV-Trolley



Flat Top Trolley

Versatile & Reliable – The Flat Top Trolley

Built for efficiency in professional settings, our Flat Top Trolley provides a durable, functional solution for transporting and organizing materials. Crafted from high-quality stainless steel with a fully welded construction, it ensures long-term reliability in demanding environments. Available in multiple sizes and shelf configurations, it features smooth swivel castors and a flat top surface for effortless loading and maneuvering, making it an essential workspace tool.

Premium Construction

Manufactured from stainless steel, this trolley is built for exceptional strength and corrosion resistance. Its fully welded structure ensures stability and durability, making it ideal for healthcare and commercial applications.

Smooth Functionality

Equipped with 4 x 100mm swivel castors, this trolley allows for effortless maneuverability, ensuring smooth movement even in confined spaces. Multiple shelf configurations provide ample storage for tools, equipment, and supplies.

Ergonomic & Practical Design

The flat top surface enables easy loading and unloading, streamlining daily tasks. Ergonomically designed to enhance workplace efficiency, it provides a stable, secure platform for transporting materials with ease.

Unmatched Quality

Designed for long-term performance, the Flat Top Trolley stands out with its robust construction and versatile applications. Whether in healthcare, hospitality, or industrial settings, its sleek, functional design ensures a reliable and professional workspace solution.



Part Number	EMY-SS20	EMY-SS20-600	EMY-SS24	EMY-SS24-600	EMY-SS26	EMY-SS28
Width	500	600	500	600	1000	500
Depth	500	500	500	500	500	500
Height	900	900	900	900	900	900
Shelves	1	1	2	2	2	3
Castors	100mm Swivel	100mm Swivel	100mm Swivel	100mm Swivel	100mm Swivel	100mm Swivel

Linen Skips

Built for stability and efficiency, this unit features precision ball-bearing castors that absorb weight and enable effortless movement. The 400mm shelf height minimizes bag lifting and helps control bag expansion, improving workplace ergonomics. Designed with a small footprint for space efficiency, it also includes full safety folds for enhanced protection. Crafted from durable, welded stainless steel, this unit is built to withstand demanding environments and is backed by a 10-Year Workmanship Guarantee, ensuring long-term reliability and industry-leading quality. 4 x 100mm swivel caster wheels.

Linen Skip - 380w x 420d x 900h

EMY-SS40LR

Linen Skip - 760w x 420d x 900h

EMY-SS42LR



IV Poles

Nesting IV Pole



Our Nesting IV Pole is designed for effortless mobility, featuring high-quality castors and height adjustability to support various medical needs. The internally weighted 5-star base enhances stability and safety, reducing the risk of tipping and pole damage. Its space-saving nesting design allows 10 poles to fit in the space of just 2, optimizing storage in busy medical facilities. Backed by a 10-Year Workmanship Guarantee, this durable, high-quality IV pole delivers long-term reliability in demanding healthcare settings.

- Size: 530w x 530d x 1400-2300h mm
- Weights: 7.2kg Internal Weights
- Castors: 5 x 75mm Swivel Castors

Nesting IV Pole EMY-SS04IWBNEST

Nylon Base IV Pole



This Nylon Base IV Pole offers height adjustability to accommodate various medical needs. Designed for stability and ease of use, it features a durable nylon base for reliable support. Built for everyday use in healthcare settings, this IV pole provides a simple, functional, and dependable solution for IV administration.

- Size: 580w x 580d x 1350-2300h mm
- Castors: 5 x 50mm Nylon

Nylon Base IV Pole EMY-SS04NB

IV Pole Accessories

Optimize your IV poles with our durable accessories, including gas bottle holders, hooks, and baskets for secure, organized storage. Designed for seamless integration, they enhance efficiency and reliability in healthcare settings.



Hooks
EMY-AO-U2



2 Hook Rams Head
EMY-RH2



4 Hook Rams Head
EMY-RH4




5 Hook Rams Head
EMY-RH5



Detachable IV Pole Hook
EMY-AO-IV



IV Pole Gas Bottle Holder



Height Adjustable Basket to Fit
on IV Poles
AO-IV-Basket 3



Height Adjustable Basket to Fit
on IV Poles
AO-IV-Basket 1/2

IV Pole Customisations

Our IV Poles can be tailored to suit various requirements, with options for welded hooks, push handles, and specialised holders for added convenience. The Nesting IV Pole can also be custom powder-coated in colors such as blue and red, offering both functionality and easy identification in busy medical environments.

Pearl White

Birch White

Lemon Yellow

Pink Blush

Flame Red

Viper Green

French Blue

Pommel Blue

Deep Violet



Push Handle
Examples



Mid Bag
Hook



Gas Bottle
Holder



Example Powder
Coat Blue



Endoscopy Travel Carts

The Featherweight® Endoscopy Cart range is designed to provide seamless transport and accessibility for essential equipment and supplies in clinical environments. With a modular design, these carts offer the flexibility to strategically position equipment and supplies for convenient access during procedures. With literally thousands of potential design configurations, each cart can be tailored to meet specific application requirements, ensuring an optimized workflow. Integrated key accessories such as scope holders, locking drawers, keyboard trays, and monitor mounts enhance functionality, while a compact footprint ensures efficient space utilization.



Featherweight® Motorized Endoscopy Cart

As a mainstay of our motorized travel carts, the Featherweight® Motorized Endoscopy Cart is designed to provide effortless transport with everything you need and nothing you don't. A retractable drive system offers flexibility, allowing users to transport and position the cart with or without power. Built for efficiency, this cart ensures smooth maneuverability and reliable performance in endoscopy suites and procedure rooms.




Featherweight® SPECS™ Motorized Endoscopy Travel Cart System

For facilities requiring a highly organized travel solution, the Featherweight® SPECS™ Endoscopy Travel Cart System introduces an innovative two-cart solution. No longer do you need to transport all your supplies on a single cart—hitch the SPECS™ trailer cart to your SPECS™ motorized equipment cart and go! This system is designed to have the smallest footprint while offering maximum organization, allowing users to separate equipment from supplies for improved workflow and accessibility during travel procedures.



Collaborative Process

Our knowledgeable representatives work with you to configure an endoscopy cart around your needed equipment. Whatever you need in an endo cart, we can make it happen.



Quality & Trust

We have been providing endoscopy teams with solutions for endoscopic equipment storage and transportation for over two decades. Your investment will last for years to come.



Configurable Design


We use modular, prefabricated components allowing each accessory, shelf, and storage option to be positioned according to your specific supply and equipment requirements.



Power Drive

Our motorized endoscopy carts are fitted with a power drive system to make travel cases a breeze. Simply engage the cart's throttle and guide the cart to its destination without strain.

Modular Design Process



Monitor Mount Options



Breakaway Monitor Cart



Monitor Arm (Options)



Monitor Stand



Monitor Riser

Accessories



Breakaway Monitor



CPU Mount



Tank Holder



Keyboard and Mouse Tray



Scope Holder



Foot Pedal Holder



Flip Up Work Surface



Wire Basket

Storage



Vertical Supply Cabinet



Tilt Bins



Drawer Dividers



Glove Box Holder



Keyless Lock



Locking Drawer Cabinet



Endoscopy Carts Designed Around You!

Endoscopy Procedure Carts

Featherweight® Procedure Room Carts keep your endoscopic equipment and supplies right where you need them to be. Our extensive range of accessories and storage systems combined with the unmatched configurability of our products means we can not only meet, but exceed, the needs of endoscopy professionals.



Double Endoscopy Cart

Featuring ample storage space, integrated cord management, and enhanced ergonomics. Configurable to your precise specifications.



Single Endoscopy Cart

A conveniently sized cart for storing supplies, imaging equipment and other endoscopic devices. An ideal solution for small procedure spaces.



Tower Endoscopy Cart

Added height allows for more storage system configurations. An excellent option for creating a small multi-role or dedicated workstation.

What We Need From You:

List of Required Equipment - eg. Light Source, Processor, Cautery
Preferred Style of Cart - eg. Procedure Room, Travel, Double, Single
Vital Accessories - eg. Scope Holder, Outlet Strips, Keyboard Tray
Essential Storage - eg. Locking Drawer, Tilt Bins, Drawer Dividers

What You Get From Us:

A custom-tailored endoscopy cart solution, designed to fit your specific equipment, storage, and workflow requirements.



Why Choose InnerSpace?

The medical industry's most effective and innovative integrated storage system



Design Configuration Options

InnerSpace specialises in designing modular storage systems (mobile carts and cabinets) bringing value to hospitals by helping organise supplies, increase efficiency and maximize storage space.

Hook Capacity

InnerSpace Heavy Duty Glide Cath Mangers, with their double track system provides maximum strength and durability. On the smooth track system which makes organising and locating catheters quickly and easy. The choice of six or eight hooks that hold can hold up to 15 catheters per hook (depending on packaging) which will maximize space and time.

Manoeuvrability

InnerSpace's heavy-duty 5 inch full swivel ball bearing plate casters (2 locking, 1 tracking) combined with a low centre of gravity provide a stable to manoeuvre cart.

Size

When space is money, InnerSpace delivers. With the best cost per cubic inch modular system. From our double wide cabinets offering an internal width of 74.9 cm and 102.9 cm for our triple wide cabinets.

Security

InnerSpace offers a range of security features on there cabinets and small carts with standard central key lock options, optional electronic keyless entry with auto re lock which takes InnerSpace level of security to another level.

Strength

In a medical environment where you really on quality from the products InnerSpace delivers with a full metal surrounding cabinet.

Price Value

Overall InnerSpace offers a total solutions build for health care. With the flexibility to change as your demands change, the choice of volume size and the overall strength, make InnerSpace cabinets modular system the only choice for medical storage.

Flexibility

There's a growing demand for storage customised to suit specialist procedure needs. It's important that storage units can be re sized and re configured to adapt to supply demands and surgical packaging changes.



An almost infinite number of combinations for tailored storage solutions

Why Choose Glass Doors?

Hinged aluminium and tempered glass doors provide 65% visibility into the cart while maintaining a low risk of supply contamination. This helps to minimise the risk of spreading infections and hinders the collection of dust within the cabinet. The added visibility also maximises efficiency, allowing staff to quickly identify and retrieve supplies by reducing the time it takes to open the cabinet.

For more information or to place an order, please contact Imaging Solutions at 1300 132 100 or visit us at www.imagingsol.com.au.

Roam Carts

Smart design. Flexible design. Durable design.

Roam's efficient exterior footprint and next-to-no wasted interior space result in more storage capacity per square inch. Roam will perform day in and day out and look good doing it.



Products shown, left to right: Roam 2 cart with accessories (SR2G) in Wilsonart Blue Agave 4919), Roam 2 counter-height cart with accessories (SR2CG).

Roam Carts are available in a wide range of configurations.



Roam 2 Arthroscopy Cart



Roam 2 Cart with 5 Shelves



Roam 2 Catheter Cart



Roam 2 Open Heart Cart



Roam 2 Ophthalmic Cart



Roam 2 Diagnostic Catheter Cart



Roam 3 Arthroscopy Cart



Roam 3 Boxed Catheter and Supply Cart



Roam 3 3-Cue Catheter and Supply Cart



See online for a full listing of
Roam Cart configurations:
[imagingol.com.au/product-
category/storage-trolleys-and-
carts/mobile-carts/roam-carts/](https://imagingol.com.au/product-category/storage-trolleys-and-carts/mobile-carts/roam-carts/)



Roam 1 counter-height cart

General storage cart configuration includes 1 single-compartment shelf, 2 shelf dividers, 1 4" basket, 1 8" basket, and label holders (SR1CGGS)



Roam 2 counter-height cart

General storage cart configuration includes 1 single-compartment shelf, 2 shelf dividers, 1 4" basket, 1 8" basket, 4 4" trays with dividers, 1 8" tray with dividers, and label holders (SR2CGGS)

Smart design

The smart design of the FlexCell system and interior components maximize Roam's storage capacity, whether you are outfitting a counter-height or full-height cart with one, two, three, or four compartments.

Glass doors are designed with a position hinge that holds the door open to exactly where you want it to be. A magnetic catch securely closes and keeps the door from opening in transit.

The roll-top door track is located on the outside of the cart, so interior storage space is not compromised.



Roam 1 cart

General storage cart configuration includes 1 single-compartment shelf, 2 shelf dividers, 2 2" trays with dividers, 3 4" trays with dividers, 2 4" baskets, 1 8" basket, 1 8" tray with dividers, and label holders (SR1RGS)



Roam 2 cart


ENT/Neuro cart configuration includes 4 single column shelves, 8 shelf dividers, 4 2" trays with dividers, 6 4" trays with dividers, 2 4" baskets, 2 8" baskets, and label holders (SR2GNU)

Durable design

Roam takes the wear and tear that mobile supply carts get. A metal base and top add stability to the case. Plastic ABS interior components won't wear down FlexCell tracks, even with frequent pulling out and pushing in of baskets and trays.

Two roll-top doors on our Roam 3 and Roam 4 make it easier to open and close and also eliminate damage common with wider roll-top doors. Both roll-top and hinged doors operate smoothly.

Even our widest Roam carts are still mobile carts, and they need to be manageable and durable. A lightweight aluminum composite frame aids maneuverability, and soft bumpers protect carts, walls, and doors against damage.



Roam 3 cart

Boxed catheter and general storage cart configuration includes 2 single-compartment shelves, 2 4" trays with dividers, 1 8" tray with dividers, 2 8" baskets, 1 double-compartment shelf, 9 shelf dividers, 1 boxed catheter module, and label holders (SR3RBG)

Flexible design

The FlexCell design and interior components create a modular system that can be configured and reconfigured as your needs and supplies change. Center columns can be retrofit to carts.

Tracks on the FlexCell system provide two tilt options for trays and baskets: a moderate tilt for most applications or a gravity-feed tilt for ER or OR applications where quick access is paramount.



Roam 4 cart

General storage cart configuration includes 8 single-compartment shelves, 16 shelf dividers, 4 2" trays with dividers, 16 4" trays with dividers, 4 4" baskets, and 4 8" baskets (SR4GGS)

Clean design

Supply carts must be highly functional, but that doesn't mean they can't look good. Concealed hinges and a flush door design add up to a better aesthetic.

Nonporous surfaces, tight seams, and few exposed fasteners support infection prevention goals and reduce dust and dirt traps.

The FlexCell system and interior components are designed for ease of cleaning. Roam carts can be cleaned with agents commonly used in hospitals. Aluminum composite surfaces won't rust or scratch.



Concealed Hinge Design

Roll-top and hinged doors can be replaced

Lightweight aluminium composite construction



Position hinge keeps door open to the position user places it



Easy-grip stainless steel push/pull handle



5" lockable medical-grade casters



**Electronic lock and keyed lock options;
no lock is standard**



**Electronic lock and keyed lock options;
no lock is standard**



Maximum Versatility

Maximize interior storage with FlexCell and a modular system of trays, baskets, and dividers.



2"-, 4"-, and 8"-high trays let you match the size of supplies to the size of the tray and reduce unused vertical space.



Stackable trays make transporting restocked trays efficient.



Powder-coated standard metal shelf spans one compartment and holds up to three cath managers®; multipurpose shelf comes in a single-compartment shelf that holds up to three cath managers and a double-compartment shelf that holds up to seven cath managers.



Labels clip on trays and baskets help organize restocking and reduce time wasted searching for supplies.



Dividers are adjustable and removable so you can organize and reorganize tray contents. Clips keep dividers securely in place.



Safety stops on the FlexCell tracks and backs of trays and baskets prevent them from sliding out of cart under normal use.



4"- and 8"-high clear baskets help you identify contents quickly. Perforated bottoms prevent dust build up.



Trays and baskets have a moderate and gravity feed tilt and are designed with straight, not sloped sides, resulting in a higher percentage of usable space and more capacity. Our 4"-high tray, compared with a sloped 4"-high tray, provides up to 38% more storage volume.

Pace Carts

Five cart heights, two widths, three drawer sizes, and a variety of accessories let you target supply storage for specific treatments and departments. A modular drawer system lets you specify the drawer configuration you need now and easily rearrange as needs and supplies change. And a host of color options lets you complement your facility's décor.



Keep pace with supply storage and transport.

Keep pace with work flow.

Keep pace with patient care.



Pace carts are mobile, easy to manoeuvre, secure, and adaptable to change.



Cart handle designed to push/pull from the side or top of handle

Antimicrobial work surface and drawer pulls

Keyed and easy-to-access, easy-to-read keypad and proximity lock options

Lightweight aluminium case

Quiet movement, easy to manoeuvre

Minimal seams, clean design



Agile and under your control.

A 17.5"-wide Pace cart is perfect for organizing supplies at the bedside or where storage space is limited. Carts with 3"- or 4"-high casters fit under counters.




Pace 21 cart with 3" casters
(SP21W4 with SCASTER3)



Prox lock (SPSDLP)

Value and performance.

Pace is designed with an antimicrobial removable surface, so you can replace this high-wear part without having to purchase a new cart. Pace carts are warranted for five years.



Antimicrobial replacement top (SPRT)

Cleanable and easy to maintain.

Made of durable aluminum, cart interiors and exteriors are designed to facilitate thorough cleaning using standard cleaners and disinfectants. Pace carts have an antimicrobial top surface and drawer pulls.



Statement of Line

**Pace 21 (SP21W4)**

Shown with 2 3/8" h, 1 6" h, and 1 9" h drawer configuration

Exterior dimensions: 29.25" w x 21.25" d x 33" h

Weight: 110 lbs

**Pace 24 (SP24W5)**

Shown with 3 3/8" h, 1 6" h, and 1 9" h drawer configuration

Exterior dimensions: 29.25" w x 21.25" d x 36" h

Weight: 120 lbs

**Pace 27 (SP21W6)**

Shown with 4 3/8" h, 1 6" h, and 1 9" h drawer configuration

Exterior dimensions: 29.25" w x 21.25" d x 39.25" h

Weight: 135 lbs

**Pace 30 (SP30W6)**

Shown with 3 3/8" h, 2 6" h, and 1 9" h drawer configuration

Exterior dimensions: 29.25" w x 21.25" d x 42.5" h

Weight: 145 lbs

**Pace 33 (SP33W7)**

Shown with 4 3/8" h, 2 6" h, and 1 9" h drawer configuration

Exterior dimensions: 29.25" w x 21.25" d x 45.5" h

Weight: 155 lbs

**Pace 21 with 3"-high casters (SP21W4 with SCASSTER3)**

Shown with 2 3/8" h, 1 6" h, and 1 9" h drawer configuration

5"-high casters are standard.

Exterior dimensions: 29.25" w x 21.25" d x 31" h

Weight: 110 lbs



Red*



Royal*



Yellow*



White*



Bisque*



Caribbean



Meadow



Lavender



Vivid Pink



Sky



Sunset

**Narrow Pace 21 (SPN21W4)**

Shown with 2 3" h, 1 6" h, and 1 9" h drawer configuration
Exterior dimensions: 17.5" w x 21.25" d x 33" h
Weight: 75 lbs

**Narrow Pace 24 (SPN24W5)**

Shown with 3 3" h, 1 6" h, and 1 9" h drawer configuration
Exterior dimensions: 17.5" w x 21.25" d x 36" h
Weight: 85 lbs

**Narrow Pace 27 (SPN27W6)**

Shown with 4 3" h, 1 6" h, and 1 9" h drawer configuration
Exterior dimensions: 17.5" w x 21.25" d x 39.25" h
Weight: 90 lbs

**Narrow Pace 30 (SPN30W6)**

Shown with 3 3" h, 2 6" h, and 1 9" h drawer configuration
Exterior dimensions: 17.5" w x 21.25" d x 42.5" h
Weight: 100 lbs

**Narrow Pace 33 (SPN33W7)**

Shown with 4 3" h, 2 6" h, and 1 9" h drawer configuration
Exterior dimensions: 17.5" w x 21.25" d x 45.5" h
Weight: 105 lbs

**Narrow Pace 21 with 3"-high casters (SPN21W4 with SCASTER3)**

Shown with 2 3" h, 1 6" h, and 1 9" h drawer configuration
Exterior dimensions: 17.5" w x 21.25" d x 31" h
Weight: 75 lbs



Red*



Royal*



Yellow*



White*



Bisque*



Caribbean



Meadow



Lavender



Vivid Pink



Sky



Sunset

Drawer sizes and combinations accommodate a variety of supplies.

Pace's modular drawer system adapts to the specific supply needs of different departments. Any combination of drawers can be specified for your cart as long as the combined height of the drawers equals the height of the cart case.



Pace 30 (SP30WF)

Shown with 4 3" (SP3DW) and 2 9" (SP9DW) drawers
Exterior dimensions: 29.25" w x 21.25" d x 42.5" h



Pace 30 (SP30WF)

Shown with 2 3" (SP3DW), 1 6" (SP6DW), and 2 9" (SP9DW) drawers
Exterior dimensions: 29.25" w x 21.25" d x 42.5" h



Pace 30 (SP30WF)

Shown with 2 6" (SP6DW) and 2 9" (SP9DW) drawers
Exterior dimensions: 29.25" w x 21.25" d x 42.5" h



3" tray with 3 long and 3 short dividers (STD3)



6" tray with 3 long and 3 short dividers (STD6)



9" tray with 1 long and 1 short divider (STD9)



3" tray with 3 long and 3 short dividers (STD3)



3" Snap-off dividers and divider clips (SSODK), includes 5 dividers and 5 clips

Interior tray dimensions Standard tray

3" tray: 21.73" w x 13.86" d x 2.47" h
6" tray: 21.73" w x 13.86" d x 5.03" h
9" tray: 21.31" w x 13.43" d x 8.15" h

Narrow tray

3" tray: 9.80" w x 13.73" d x 2.47" h
6" tray: 9.66" w x 13.60" d x 5.03" h
9" tray: 9.50" w x 13.43" d x 8.15" h

Organize drawer interiors with trays.

Three tray heights are sized to drawer heights to maximize every inch of space. Trays lock into place to keep contents from shifting when carts are wheeled over door and elevator thresholds.

Further organize supplies with adjustable tray dividers.



Trays are sized to fit both standard width and narrow Pace carts.

Shown with 6"h narrow tray with 1 long and 3 short dividers (SNTD6)



Optional tray lid (SSTL6) for 6"h smooth tray (SST6)



Pace transport carts have FlexCell interiors and open fronts. Trays, drawers, and baskets slide into the cells to transport supplies within a department or among floors.



Products shown, from left: Pace transport cart SP27WC with 1 3"h (ST3), 1 6"h (ST6), and 2 9"h trays (ST9) and label holders (SLH35),

Pace 30 cart with Pharmacy configuration (SP30W6 with SPPVP)

Accessories define the function of your cart.

From glove box holders and waste containers on carts used for general patient care to monitor arms to hold displays on carts in an imaging department, accessories increase the functionality of every cart.

Attaching and removing accessories is easy and quick, so changing or adding to the functionality of your cart is a snap.



Products shown: Pace 27 cart with Anesthesia configuration (SP27B6 with SPAVP)



A slide-out surface extends from either side of a standard Pace cart to create work and storage space; the slide-out surface for a narrow Pace cart slides to the left side.

Shown with optional narrow pull-out work surface (SPNPOWS)



Shown with optional pull-out work surface (SPPOWS)

Lock options keep supplies secure.



Keypad Lock (SPSDLE)



Proximity Lock (SPSDLP)



Pace carts are offered with several lock options, including keypad, badge swipe, and keypad and badge swipe. Keyed lock is standard. Proximity locks can be programmed and updated via computer and USB connection.

Products shown: Pace 27 cart with Code configuration (SP27R6 with SPCVP)



Breakaway lock bar (SPLB)



Single drawer breakaway lock (SPSBL)

Specialty carts. Ready with a purpose in mind.

Code carts, isolation carts, anesthesia carts, difficult airway carts, cast carts, IV prep carts, pharmacy carts, pediatric care carts, and transport carts come preconfigured and accessorized. Specialty carts can be specified in any height of a standard 29.25"-wide Pace cart.* Transport carts are offered on both standard and narrow carts.

*Difficult airway cart not available for 21"-high cart



Pace Pharmacy Cart

(SP21W4, SP24W5, SP27W6, SP30W6, and SP33W7 with SPPVP)

Pace 21

Exterior dimensions: 37.25"w x 21.25"d x 39"h

Pace 24

Exterior dimensions: 37.25"w x 21.25"d x 42"h

Pace 27

Exterior dimensions: 37.25"w x 21.25"d x 45.25"h

Pace 30

Exterior dimensions: 37.25"w x 21.25"d x 48.25"h

Pace 33

Exterior dimensions: 37.25"w x 21.25"d x 51.5"h

Pace Anaesthesia Cart

(SP21B4, SP24B5, SP27B6, SP30B6, and SP33B7 with SPAVP)

Pace 21

Exterior dimensions: 42"w x 21.6"d x 57.25"h

Pace 24

Exterior dimensions: 42"w x 21.6"d x 60.5"h

Pace 27

Exterior dimensions: 42"w x 21.6"d x 63.5"h

Pace 30

Exterior dimensions: 42"w x 21.6"d x 66.75"h

Pace 33

Exterior dimensions: 42"w x 21.6"d x 69.75"h



Pace IV Prep Cart

(SP21W4, SP24W5, SP27W6, SP30W6, and SP33W7 with SPIVP)

Pace 21

Exterior dimensions: 37.25"w x 21.6"d x 59.75"h

Pace 24

Exterior dimensions: 37.25"w x 21.6"d x 63"h

Pace 27

Exterior dimensions: 37.25"w x 21.6"d x 66"h

Pace 30

Exterior dimensions: 37.25"w x 21.6"d x 69.25"h

Pace 33

Exterior dimensions: 37.25"w x 21.6"d x 72.5"h



Pace Code Cart

(SP21R4, SP24R5, SP27R6, SP30R6, and SP33R7 with SPCVP)

Pace 21

Exterior dimensions: 39.6"w x 21.6"d x 59.75"h

Pace 24

Exterior dimensions: 39.6"w x 21.6"d x 63"h

Pace 27

Exterior dimensions: 39.6"w x 21.6"d x 66"h

Pace 30

Exterior dimensions: 39.6"w x 21.6"d x 69.25"h

Pace 33

Exterior dimensions: 39.6"w x 21.6"d x 72.5"h



Pace Isolation Cart

(SP21Y4, SP24Y5, SP27Y6, SP30Y6, and SP33Y7 with SPISVP)

Pace 21

Exterior dimensions: 34.6"w x 21.25"d x 33"h

Pace 24

Exterior dimensions: 34.6"w x 21.25"d x 36.25"h

Pace 27

Exterior dimensions: 34.6"w x 21.25"d x 39.25"h

Pace 30

Exterior dimensions: 34.6"w x 21.25"d x 42.5"h

Pace 33

Exterior dimensions: 34.6"w x 21.25"d x 45.5"h



Pace Cast Cart

(SP21W4, SP24W5, SP27W6, SP30W6, and SP33W7 with SPCTVP)

Pace 21

Exterior dimensions: 37"w x 21.6"d x 57.25"h

Pace 24

Exterior dimensions: 37"w x 21.6"d x 60.5"h

Pace 27

Exterior dimensions: 37"w x 21.6"d x 63.5"h

Pace 30

Exterior dimensions: 37"w x 21.6"d x 66.75"h

Pace 33

Exterior dimensions: 37"w x 21.6"d x 69.75"h



Pace Paediatric Care Cart

(SP27PED)

Exterior dimensions: 29.25"w x 21.25"d x 39.25"h



Shown with 9 optional individual breakaway locks (SPSBL)



Pace Difficult Airway Cart

(SP24W5, SP27W6, SP30W6, and SP33W7 with SPDVP)

Not available on a Pace 21 cart

Pace 24

Exterior dimensions: 50.6"w x 21.6"d x 63"h

Pace 27

Exterior dimensions: 50.6"w x 21.6"d x 66"h

Pace 30

Exterior dimensions: 50.6"w x 21.6"d x 69.25"h

Pace 33

Exterior dimensions: 50.6"w x 21.6"d x 72.5"h



Pace Isolation Cart

(SP21Y4, SP24Y5, SP27Y6, SP30Y6, and SP33Y7 with SPISVP)

Pace 21

Exterior dimensions: 34.6"w x 21.25"d x 33"h

Pace 24

Exterior dimensions: 34.6"w x 21.25"d x 36.25"h

Pace 27

Exterior dimensions: 34.6"w x 21.25"d x 39.25"h

Pace 30

Exterior dimensions: 34.6"w x 21.25"d x 42.5"h

Pace 33

Exterior dimensions: 34.6"w x 21.25"d x 45.5"h



Pace Cast Cart

(SP21W4, SP24W5, SP27W6, SP30W6, and SP33W7 with SPCTVP)

Pace 21

Exterior dimensions: 37"w x 21.6"d x 57.25"h

Pace 24

Exterior dimensions: 37"w x 21.6"d x 60.5"h

Pace 27

Exterior dimensions: 37"w x 21.6"d x 63.5"h

Pace 30

Exterior dimensions: 37"w x 21.6"d x 66.75"h

Pace 33

Exterior dimensions: 37"w x 21.6"d x 69.75"h

Turn your Pace cart into a workstation on wheels with a height-adjustable surface accessory.

- 12" pneumatic height-adjustment range
- Platform holds a laptop or standard-sized keyboard
- Raised-edge surfaces
- Available for all heights of standard-width and narrow carts
- Adjustable work surface can be specified preconfigured on Pace 21 standard-width and narrow carts



A sampling of the many accessories available for Pace carts.



Chart Holder (SCHP)



Glove Box Holder, single (S1GBHP)



Wire Basket (SSMWBP)



Multi-Roll Tape Dispenser (SSTD)



Flip-Up Work Surface (SFWSP)



Sharps Container (SSCP)



Pull-Out Keyboard Tray (SPKBT)



Push Handle (SPAPH)



Security Lock Box, narrow (SPNSLB3)



Tilt Bin, 4 compartment (STB4)



Top-Mounted Accessory Rail, double (SPUAR2)



Pegboard Hooks (SRPB4H, SRPB6H, SRPB8H, SRPB11H)



Rotating SneezeGuard Shield (SPSGS)



Adjustable Work Surface, narrow (SPNAWS)



Side-Mounted Scope Cabinet (SSMSCP)



IV Pole (SIVPP)



Surface Rail (SPWAR)




Defibrillator Shelf (SDSP)



Flat Panel Monitor Arm (SFPMAP)



Disinfectant Holder (SRDHP)



Catheter Holder (SPCH)



Back-Mounted Adjustable Rail (SAMRBP)



Cord Wrap (SCWP), Power Strip (SPSP)



Resuscitation Board (SRBP)

RADstore Supply Carts

The medical industry's most effective and innovative integrated storage system, now assembled in Australia!

Designed for versatility, durability, and affordability, RADstore by Imaging Solutions is an Australian-made, cost-effective alternative to traditional supply carts. Developed in collaboration with Solaire, RADstore maintains full compatibility with the market-leading FlexCell modular storage system, ensuring seamless integration with a wide range of storage solutions.

Smart Design

RADstore carts maximize storage efficiency with the Solaire FlexCell system, offering an optimized interior layout that enhances capacity while ensuring easy access to supplies. This space-saving design makes it an ideal solution for healthcare, radiology, and general storage applications.

Durable Design

Built for daily use, RADstore is constructed from high-quality laminated MDF with a reinforced base for added stability. Its ABS interior components provide smooth, wear-resistant operation, ensuring FlexCell trays and baskets remain reliable, even with frequent restocking and adjustments.

Flexible Design

The FlexCell modular system enables customized storage configurations that adapt as your needs evolve. With a comprehensive range of accessories RADstore with FlexCell, offers a scalable and organized storage solution the medical environment.

Clean Design

RADstore carts feature a seamless, nonporous finish that supports infection control by minimizing dust and dirt traps. The laminated surfaces and FlexCell components are easy to clean with standard hospital-grade disinfectants, ensuring long-term hygiene and durability in clinical settings.

Example Configurations



Mixed Catheter Cart



Diagnostic Catheter Cart



Boxed Catheter Cart



Suture Cart

Customised Colours, Decals and Cobranding Available On Request



Narrow Cart

RADstore's Narrow Cart offers a compact yet highly efficient storage solution, ideal for tight spaces while maintaining full compatibility with Solaire's FlexCell modular system. Designed for durability and ease of access, it provides a customizable storage layout to suit a variety of applications.

- Height: 42 FlexCells
- Width: 1 Compartment
- Dimensions: 620 x 730 x 2050mm (LxDxH)

Narrow Cart, Single Compartment

RAD-ST-SC1R

Double Cart

RADstore Double Cart delivers maximum storage capacity with the option of two compartments for enhanced organization or a single wide compartment for flexible configurations. Built for durability and efficiency, it seamlessly integrates with Solaire's FlexCell trays, baskets, and accessories, offering a versatile, high-performance storage solution.

- Height: 42 FlexCells
- Width: 2 Compartments or 1 Wide Compartment
- Dimensions: 1058 x 730 x 2050mm (LxDxH)

Double Cart, Dual Compartments

RAD-ST-SC2R

Double Cart, Wide Compartment

RAD-ST-SC2RNC

Accessories



Shelves

Shelves for easy access storage, optional dividers with notches designed to hold backstops for smaller boxes. Shelves can also be used to hang cath managers.

Plastic Shelf

- Height: 1 FlexCell, 4 with Divider
- Width: 1 Compartment
- Storage: 23.5" w x 16" d
- Weight: Up to 11 kg (25 lbs)

Plastic Shelf SLM-S27P1

Divider SLM-S27DI

Divider Backstop SLM-S27BS1

Metal Shelf

- Height: 1 FlexCell, 4 with Divider
- Width: 1 Wide Compartment
- Storage: 23.5" w x 16" d
- Weight: Up to 22 kg (50 lbs)

Metal Shelf SLM-S27M2

Divider SLM-S27DI

Divider Backstop SLM-S27BS2



Cath Manager Module

Module includes a shelf and cath managers with each inline hook holds 10-15 catheter packages. Cath managers can be easily adjusted horizontally, add or removed.

1 Compartment, Cath Manager Module

- Height: 5 FlexCell, minimum
- Width: 1 Compartment
- Storage: 3 x 9-Hook on Plastic Shelf

3x 9-Hook Module SLM-S27P1CM3

2 Compartment, Cath Manager Module

- Height: 5 FlexCell, minimum
- Width: 1 Wide Compartment
- Storage: 7 x 9-Hook on Metal Shelf

7x 9-Hook Module SLM-S27M2CM7

Accessories / Individual

1x 9-hook Cath Manager SLM-S9CM

Label Cards (25 Pack) SLM-S9CMCARDS



Trays

Tray has seamless construction for ease of cleaning. Features an integrated safety backstop while maintaining accessibility for cleaning or reconfiguration.

6"H Tray, Grey

- Height: 4 FlexCell
- Width: 1 Compartment
- Storage: 13.86" w x 21.73" d x 5.03" h
- Weight: Up to 18 kg (40 lbs)

6"H Tray SLM-ST6

6"H Clear Short Adj. Divider SLM-SSD6

6"H Clear Long Adj. Divider SLM-SLD6

9"H Tray, Grey

- Height: 6 FlexCell
- Width: 1 Compartment
- Storage: 13.86" w x 21.73" d x 5.03" h
- Weight: Up to 18 kg (40 lbs)

9"H Tray SLM-ST9

9"H Clear Short Adj. Divider SLM-SSD9

9"H Clear Long Adj. Divider SLM-SLD9



Boxed Catheter Insert

Powder-coated steel box and base with two pull-out sections for double wide carts. Each section contains two adjustable clear acrylic dividers for a total of 6 inserts. Pull-out sections feature anti-tip mechanism for added safety and usability.

- Height: 29 FlexCells
- Width: 1 Wide Compartment
- Storage: 15" w x 21" d x 38.24" h.
- Weight: Up to 22 kg (50 lbs)
- Extends: approx. 12.75" from edge

Boxed Catheter Insert SLM-S27BC2-AT



Basket, Clear

Perforated construction keep dust out of baskets. Features an integrated safety backstop while maintaining accessibility for cleaning or reconfiguration.

6"H Basket, Clear

- Height: 4 FlexCell
- Width: 1 Compartment
- Storage: 13.86" w x 21.73" d x 5.03" h
- Weight: Up to 18 kg (40 lbs)

6"H Basket, Clear SLM-SB6

6"H Clear Short Adj. Divider SLM-SSD6

6"H Clear Long Adj. Divider SLM-SLD6

9"H Basket, Clear

- Height: 6 FlexCell
- Width: 1 Compartment
- Storage: 13.43" w x 21.31" d x 8.15" h
- Weight: Up to 18 kg (40 lbs)

9"H Basket, Clear SLM-SB9

9"H Clear Short Adj. Divider SLM-SSD9

9"H Clear Long Adj. Divider SLM-SLD9



Suture Module

Module includes 6 pull-out shelves to store and access suture boxes and other small supplies. Each shelf holds 8x 2.25"-deep suture boxes for a total of 48 boxes. Module is expandable to 12 shelves (6 per side)

- Height: 18 FlexCells
- Width: 1 Compartments

Suture Module SLM-SPSM

Suture Shelf Dividers SLM-SSSHDI

Additional Shelves SLM-SSSH

Label Holders

Label holder for 2" x 3.5" cards clips to the front of all trays and baskets. Label holder for 1.5" x 2" cards clips to an interior tray or basket divider. Cards not included.

Dividers (Pack 10) SLM-SLH25

Trays/Baskets (Pack 10) SLM-SLH35

Explore Our Range


At Imaging Solutions, we take pride in offering an extensive and comprehensive range of products designed to cater to the diverse needs of healthcare facilities across the globe. Our commitment to providing top-quality equipment and technology solutions is at the core of our Preferred Supplier Agreement. This ensures that our customers receive unparalleled access to the best in class products, while benefiting from the tangible value we deliver through managing and containing operational costs. With a broad spectrum of offerings, our Single Source Supply agreements are the ideal solution for meeting your facility's specific needs.

As your trusted partner, Imaging Solutions is dedicated to providing a seamless and integrated experience, allowing you to focus on delivering exceptional patient care. With our Preferred Supplier Agreement, you can have confidence that you are investing in the most advanced and reliable products in the industry, all while enjoying the cost-saving benefits that our Single Source Supply agreements provide. Explore Our Range and discover the Imaging Solutions difference.


Your Single Source Supplier™

Established initially as a small medical imaging hardcopy film and consumables distributor, specialising in analogue based imaging technology, the company very quickly identified the need for a reliably supply source of high quality imaging accessories. As a result the business targeted a small niche of premium suppliers of imaging accessories and over a short space of time acquired distribution access these to develop an initial foundation for today's expansive and comprehensive accessories portfolio. These industry wide market leading brands converge together to develop a comprehensive product portfolio offering the company markets today.


Effectively on completion of our initial growth phase, in building the accessories range offering, Imaging Solution set about the task of identifying the world's market leading brand, in every functional specialisation existing in the imaging market.



Radiation Protection



Positioning and Support Aids




Grids and Detector Accessories




MRI Equipment and Accessories




Design and Construction




Patient Experience



Imaging Markers and Procedure Accessories




Gels, Contrast and Warmers




Tables, Stools, Transport and Furniture



Trolleys, Stands, Carts and Storage



Suspension Systems




Surgical and Examination Lighting




Display Systems



Data Management



Healthcare IT



Phantoms

your single source supplier™

your single source supplier™

Global / Online

www.imgsol.com
info@imgsol.com
+61 7 3387 0400

Australia

PO Box 3225
Loganholme, QLD, 4129
Australia
1300 132 100

New Zealand

PO Box 113098
Newmarket, Auckland
New Zealand
0800 723 776

United States

222 W. Las Colinas Blvd, Suite 1650E
Irving, TX, 75039
USA
(866) 586-0915

Europe

25 Bd Romain Rolland,
92120 Montrouge
France
+33 6 5262 3244

